

lifetime limited warranty

Your Drive branded product is warranted to be free of defects in materials and workmanship for the lifetime of the product for the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.



© 2015 Medical Depot, Inc. All rights reserved.
Drive is a trademark of Medical Depot, Inc.
Port Washington N.Y. 11050 USA Made in China

Rev1.1.09.15

cane seat

item # 10365



EU Authorized Representative
Drive Medical LTD
Ainley's Industrial Estate
Elland, West Yorkshire,
United Kingdom HX5 9JP

www.drivemedical.com

A large, vertical version of the Drive logo, with the word "drive" in white lowercase letters on a dark blue background. A small yellow and orange square is positioned above the letter 'i'. A trademark symbol (TM) is located at the top right of the word.

cane seat



Item # 10365 Shown



WEIGHT CAPACITY: 250 LBS (113 KG)

The Drive Cane Seat is a versatile and convenient seating device, when opened and a strong support cane when closed. This device is intended for individuals who need to rest when walking a long distance. It is equipped with a durable plastic seat. Individuals can comfortably hold onto this device on the padded foam grip handles.

features

- No assembly required
- Unfolds easily in seconds
- Durable aluminum bronze construction
- Cushioned foam grip handles
- Comfortable plastic seat

proper sitting procedure



lifetime limited warranty

Your Drive branded product is warranted to be free of defects in materials and workmanship for the lifetime of the product for the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.



© 2014 Medical Depot, Inc. All rights reserved.

Drive is a trademark of Medical Depot, Inc.

Port Washington N.Y. 11050 USA Made in China

Rev.5.05.20.14

adjustable height cane seat

item # RTL10365-ADJ



EU Authorized Representative
Drive Medical LTD
Ainley's Industrial Estate
Elland, West Yorkshire,
United Kingdom HX5 9JP

www.drivemedical.com



cane seat



Item # RTL10365-ADJ Shown

WEIGHT CAPACITY: 250 LBS (113 KG)

The Drive Cane Seat is a versatile and convenient seating device when opened and a strong support cane when closed. This device is intended for individuals who need to rest when walking a long distance. It is equipped with a durable plastic seat. Individuals can comfortably hold onto this device on the contoured grip handles.

features

- No assembly required
- Unfolds easily in seconds
- Durable aluminum bronze construction
- Contoured foam grip handle
- Comfortable plastic seat

proper sitting procedure

