



SAFEHIP® AirX™

- the 100% textile hip protector



Excellent user comfort

- The textile protectors are hardly visible
- Flexible and comfortable pants with a perfect fit
- Can be used with pads for light/medium incontinence
- Pants have no seams and no labels

Excellent breathability

- Minimises skin irritations such as reddening and itching
- Gives a dry and comfortable feeling

Washable

- 95°C – pants with protectors
- Fast drying, tumble dryer safe

SAFEHIP® size tables

Measure the dimensions of the broadest part of your hips to find out the appropriate size. If the dimensions are between two sizes, the smaller size should be chosen if the user has slim thighs. The larger size should be chosen in the case of more sturdy thighs.

SAFEHIP® AirX™

SAFEHIP® Pants	Cm	In	Measure your hips
XS	65-85	26-34	
S	75-95	30-38	
M	90-110	36-44	
L	100-120	40-48	
XL	110-140	44-56	
XXL	120-150	48-60	

SAFEHIP® Active

SAFEHIP® belt	Cm	In	Measure your hips
S	85-95	34-38	
M	95-105	38-42	
L	105-115	42-46	
XL	115-125	46-50	
XXL	125-140	50-56	



B-Safehip-Enduser-Leaflet-EN-01-0515

CE www.SAFEHIP.com

Tytex A/S · Industrivej 21 · 7470 Ikast
mail@tytex.com · Tel. +45 9660 4200



**Break the fall,
not the hip!**

