

1. Knee Pad Screws (4 ea)- 95012R79601
2. Brake Cable Only- 95012R79602
3. Wheel Lock w/o Pad Left- 95012R79603
Right- 95012R79604
4. Wheel Lock Bolts- 95012R79627
5. Wheel Lock Spring- 95012R79605
6. Brake Connecting Bar- 95012R79606
7. Brake Pad & Hardware- 95012R79615
8. Front or Rear Axle- 95012R79608
9. Front or Rear Wheel Assm.- 95012R79609
10. Hand Brake & Cable- 95012R79610
11. Knee Pad Post- 95012R79611
12. Knee Pad & Hdwr (ea 2 needed)- 95012R79612
13. Seat Post Release Lever- 95012R79613
14. Knee Post Adjustment Pin- 95012R79614
15. Tiller Adjustment Knob- 95012R79607
16. Fork Stem Bearing- 95012R79616
17. Hand Grip (ea)- 95012R79617
18. Cable Retaining Ring- 95012R79618
19. Steering Handle Folding Hinge- 95012R79619
20. Steering Column Tension Spring- 95012R79620
21. Handle Bar- 95012R79621
22. Front Fork- 95012R79622
23. Handle Bar Basket Mount- 95012R79623
24. Folding Release Lever- 95012R79624
25. Axle Retaining Bolts (4 ea)- 95012R79625
26. Basket- 95012R79626
27. Wheel Bearing (ea)- 9503H79021

