

Steerable Knee Walker



- Knee walker can be steered for increased maneuverability
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot
- New, two-piece padded leg pad provides comfort
- Deluxe braking system
- Removable front basket
- Durable, steel frame with silver vein finish
- Tool free height adjustment
- 7" casters are ideal for indoor/outdoor use

PRODUCTS**Description**

796	Economy Knee Walker, 1/cs
-----	---------------------------

SPECIFICATIONS**Description**

Handle Height	31"-40"
Unit Dimensions	16"-30"
Knee Pad Height	16.5"-20.5"
Pad Dimensions	7" (W) x 15" (D)
Weight	22 lbs.
Carton Shipping Weight	27 lbs.
Weight Limit	350 lbs.
Warranty	Limited Lifetime, Lifetime Warranty on Brake Cable